

# Midmark 626 Power Exam Table - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Midmark 626 Power Exam Table - Refurbished. The Midmark 626 Barrier-Free Examination Chair enhances accessibility with a low seat height and flat transfer surface. Patients can easily transfer on and off the chair with minimal help. Its design improves efficiency and promotes accurate blood pressure readings during exams.

Clinics benefit from powered back movement, pelvic tilt, and quick exam functions that save time and improve comfort. This Midmark 626 exam chair includes hand and foot controls, stirrups, accessory receivers, and a heated drawer. Clean Assist Roller System allows easy repositioning for sanitation.

Tamper-resistant receptacles on both sides support safe power use for medical devices.

## FEATURES

- Easy patient positioning for better blood pressure measurements.
- Powered movement of the back section ensures the patient's back is supported.
- Barrier-Free low chair height allows the majority of patients, even female patients with heights in the third percentile, to place their feet flat on the floor.
- The Midmark IQconnect framework provides a seamless connection from your devices to your EMR system. In other words, it allows you to operate all of your devices and transfer data to the EMR from one single "pane of glass" -- your computer screen.
- Clean Assist Roller System allows users to easily and safely move the chair for cleaning.

## SPECIFICATIONS

- Patient Weight Capacity: 650 lb (295 kg).
- Height
- Uncompressed: 15.5" Low / 37" High (39 cm / 94 cm).
- Compressed: 14.6" Low (36 cm).
- Length with Footrest Extended: 76" (193 cm).
- Pelvic Tilt: 5.2°.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)