

Midmark 204 Exam Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Midmark 204 is a manual exam table designed for patient comfort and clinical efficiency. It features a manual backrest with air spring assist for smooth positioning, an enlarged foot section 40% larger than standard models for enhanced patient stability, and three pass-through drawers accessible from either side. The one-piece seamless design prevents spills and supports quick cleaning to maintain sanitary conditions.

The table accommodates a 500 lb patient weight capacity and includes a paper roll holder for standard 21" x 3.5" rolls. Optional accessories include a 7° pelvic tilt and 35-watt heated drawer for patient warmth. A duplex hospital-grade electrical receptacle provides power convenience for medical devices.

Constructed with powder-coated steel and smooth drawer slides, the Midmark 204 is built for durability. It comes completely refurbished with new upholstery available in multiple colors and finishes. The table meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 electrical safety standards.

FEATURES

- 500 lbs. patient weight capacity for universal use
- New upholstery in your choice of colors
- Seamless drawers with rounded edges for infection control
- Enlarged foot section for safer patient access
- Manual back with air spring assist for smooth positioning
- Durable steel frame with powder-coated finish
- Easy-to-clean antimicrobial upholstery
- Optional accessories for specialty procedures

SPECIFICATIONS

- Maximum Patient Weight: 500 lb (226.8 kg).
- Weight of Table: 220 lb (99.8 kg).
- Height with Upholstery: 33.13" (84.2 cm).
- Seamless Upholstery Top: 28.5" (W) x 59.5" (L) (72.4 cm x 151.1 cm).
- Paper Roll Holder: 21" x 3.5" (53.3 cm x 8.9 cm).
- Length with Footrest Extended: 76" (193 cm).
- Electrical Receptacle: 115 VAC, 5 AMPS Maximum; Duplex, hospital grade electrical receptacle on left side.
- Pelvic Tilt: 1.5° or 7° (optional).

**Brochure and images are a general overview of the product. Actual product options may vary.*