

# Midmark Ritter 75E Evolution Power Procedure Chair/Table - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Midmark Ritter 75E Evolution Power Procedure Chair/Table - Refurbished. The refurbished Midmark Ritter 75E Evolution Power Procedure Chair/Table offers a dependable and cost-effective solution for today's medical practices.

For decades, healthcare professionals have trusted Midmark for high-quality, durable equipment that enhances patient care and operational efficiency. The 75E Evolution continues that tradition with robust features and reliable performance in a wide range of clinical settings.

The Midmark Ritter 75E is ideal for general exams, dermatology, OB/GYN, and outpatient procedures. The Midmark 75E has a footswitch that lets you control all powered movements without using your hands.

## FEATURES

- Movements: Table Up/Down, Tilt Foot Up/Down, and Back Tilt Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper Holder

## SPECIFICATIONS

- Minimum Height - 26" (66 cm)
- Maximum Height - 42" (106.7 cm)
- Width (top) - 27" (68.6 cm)
- Weight - 475 lbs.
- Tilt - Trendelenburg to 45 Degrees

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)