


# Mindray Passport 8 Patient Monitor with CO2 - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Mindray Passport 8 is a compact patient monitor featuring an 8-inch touchscreen display with user-friendly touch function keys. It provides essential monitoring parameters including 3/5-lead ECG, SpO2, NIBP, respiration, temperature, and CO2 measurement.

The monitor's modular plug-and-play design allows for easy expansion with advanced monitoring options such as invasive pressure and cardiac output. It supports 96 hours of graphic and list trends, visual and audible alarms, and remote viewing with multi-bed alarm notifications.

Device connectivity is available via DIAP and Ethernet, and the unit is powered by a rechargeable lithium-ion battery. The package includes a CO2 module, new ECG patient cable and lead wires, new SpO2 finger sensor, new NIBP cuff and hose, new battery, power cord, printer paper, and user manual.

## FEATURES

- 8.8" L x 5.3" W x 9.5" H Chassis

## SPECIFICATIONS

- 8.8" L x 5.3" W x 9.5" H Chassis

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)