

Mindray Passport 12 Patient Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Mindray Passport 12 Patient Monitor – Refurbished. The Mindray Passport 12 is a versatile, modular patient monitor designed to adapt seamlessly to various clinical environments. Additionally, its plug-and-play modular architecture allows for advanced monitoring, including CO2, dual invasive pressures, and cardiac output.

Furthermore, an optional multi-gas module automatically identifies and measures gas concentrations. Also, this includes O2, CO2, N2O, and five anesthetic agents. Moreover, this monitor offers integrated HL7 and DIAP communication protocols, with both wired and wireless networking options.

As a result, it connects to the BeneVision Central Station and Panorama Patient Monitoring System for centralized patient data management.

FEATURES

- 12" high-resolution touchscreen display
- Up to 8 waveforms
- Standard features include 3 or 5-lead ECG, Masimo SET® SpO2, 2 IBP, NIBP, respiration, temperature, integrated recorder
- Optional Nellcor® OxiMax® SpO2 available
- Available ST and arrhythmia analysis

SPECIFICATIONS

- Monitor Size: 10.9 H x 5.4 D x 12.7 W in
- Weight: 10 lb (standard parameters configuration, including 2 lithium batteries and a recorder)
- Power Supply: AC power, battery
- AC Power: 100 to 240 VAC (±10%)
- Current: 1.3 to 0.5 A

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)