

Mindray Passport 12m Patient Monitor - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Mindray Passport 12m is a patient monitor designed for intermediate to critical care settings. It features a 12.1-inch high-resolution touchscreen display with an intuitive interface and customizable screen views that allow healthcare professionals to quickly access relevant patient data. The monitor supports quick-access menus and touch function keys for efficient operation.

The Passport 12m provides comprehensive monitoring capabilities including multi-lead arrhythmia analysis, 12-lead ECG interpretation, and support for up to 8 invasive blood pressure measurements. Optional modules such as multi-gas and BIS technologies can be added to expand functionality based on clinical requirements.

The monitor integrates seamlessly with EMR systems for accurate data export and streamlined workflows in both inpatient and operating room settings. It offers extensive data storage capability for trend data, alarms, events, and 48 hours of full disclosure recording.

FEATURES

- 12.1" high-resolution touchscreen display.
- Familiar intuitive interface across the Passport family.
- Customizable screen views and quick keys for diverse patient and clinical settings.
- Robust measurement set for intermediate to critical care environments.
- Supports T1 Transport Monitor/Module as a multi-parameter module.
- Extensive data storage capability for trend data, alarms, events, and 48 hours of full disclosure.

SPECIFICATIONS

- Height: 11.7 in (29.7 cm)
- Width: 13.2 in (33.6 cm)
- Depth: 7.4 in (18.7 cm)
- Weight: 16.5 lbs (7.5 kg) (Without modules, batteries, accessories, and recorder)
- Size (Diagonal): 12 in (30.5 cm)
- Resolution: 800 x 600
- Run Time (Average): About 330 minutes when powered by two new fully-charged batteries
- Charge Time: Less than 6 hours to 100%

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)