

Mindray Passport 12m Patient Monitor w/ AG Module - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Mindray Passport 12m Patient Monitor w/ AG Module - Refurbished. The refurbished Mindray Passport 12m combines clinical accuracy with flexibility, making it ideal for mid to high-acuity care settings. Featuring an integrated AG (Anesthesia Gas) module, it provides continuous monitoring of anesthesia gases for enhanced respiratory and anesthetic management.

Its 12.1-inch touchscreen display offers clear visualization and intuitive control. Quick-access menus and touch keys simplify operation, helping clinicians maintain focus on patient care during critical situations. This monitor includes important features like arrhythmia detection and 12-lead ECG monitoring.

It can track up to eight blood pressure lines and connects easily to EMR systems.

FEATURES

- 12.1" high-resolution touchscreen display.
- Familiar intuitive interface across the Passport family.
- Customizable screen views and quick keys for diverse patient and clinical settings.
- Robust measurement set for intermediate to critical care environments.
- Supports T1 Transport Monitor/Module as a multi-parameter module.

SPECIFICATIONS

- Height: 11.7 in (29.7 cm)
- Width: 13.2 in (33.6 cm)
- Depth: 7.4 in (18.7 cm)
- Weight: 16.5 lbs (7.5 kg) (Without modules, batteries, accessories, and recorder)
- Type: Color TFT Display

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)