

# MPD Collapsible Side Shelf Cart Accessory | TSN-2



*Stock Photo - Options May Vary*

## **FACTORY NEW**

MPD Cart Accessory: Collapsible Side Shelf TSN-2. Adds valuable workspace to your Standard Series medical cart. Designed for flexibility, this shelf offers extra surface area when needed and folds down easily when not in use.

Its compact size—17-3/4" deep by 14-7/8" wide—makes it a practical solution for tight clinical spaces. MPD's collapsible side shelf supports a wide range of tasks, from documentation to temporary supply placement.

It installs securely on Standard Series Carts, offering dependable support without taking up drawer space. With its space-saving design, the TSN-2 is ideal for environments where mobility and efficiency are essential.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)