

# MPD Cart Defib Swivel Shelf Unit Accessory | TDS-1



*Stock Photo - Options May Vary*

## FACTORY NEW

MPD Cart Accessory: Defib Swivel Shelf Unit TDS-1. The MPD TDS-1 adds versatility and function to your medical cart setup. In addition, this swivel defib shelf measures 14-1/16" deep by 26-7/8" long, providing reliable support for defibrillator units.

The top shelf rotates a full 360°, allowing quick, efficient access from any angle during critical procedures. Additionally, this single swivel shelf adds 12" of height to your MPD cart, optimizing space without sacrificing mobility or convenience. Furthermore, it installs easily and offers a secure platform for essential emergency equipment.

For best results, pair it with a compatible MPD cart. Please note: Cart not included, and is sold separately.

## FEATURES

- Swivel Defib Single Shelf
- Adds additional 12" to the height of the cart

## SPECIFICATIONS

- Model Number: TDS-1 Swivel defibrillator single shelf measures:
- Dimensions: 14-1/16" Deep x 26-7/8" Long
- Swivel Range: Top swivels 360°

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)