

MPD Mini/Tower Full Drawer Tray Cart Accessory | TMH-10



Stock Photo - Options May Vary

FACTORY NEW

The TMH-10 Full Drawer Tray is a versatile accessory designed for MPD Mini and Tower Carts to simplify drawer organization and improve medication storage. It features 4 dividers, 10 ampule holders, and 2 rails that work together to keep contents organized and secure during cart movement.

This tray helps reduce clutter while keeping supplies neat, visible, and readily accessible. The secure design ensures items remain in place even with frequent cart repositioning, supporting efficient workflow in clinical, medical, and laboratory environments.

The TMH-10 enhances daily operations by keeping essential items within reach and organized. It can be customized to match specific storage needs and integrates seamlessly with MPD Mini and Tower Cart systems.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)