

MPD Optional Adjustable Fluid Tray System Cart Accessory | TFT-24



Stock Photo - Options May Vary

FACTORY NEW

MPD Cart Accessory: Optional Adjustable Fluid Tray System TFT-24. Upgrade your cart's organization and storage efficiency with the TFT-24 Optional Adjustable Fluid Tray System by MPD. Designed specifically for 24" space carts, this system provides a practical solution for neatly storing fluids, bottles, and other medical supplies.

Moreover, it's ideal for healthcare professionals who need customizable, high-capacity storage in fast-paced environments. Each TFT-24 package includes three adjustable fluid trays, giving you flexibility to arrange items as needed.

Additionally, the TFT-24 tray, shown on the far left, measures 18-3/8" wide x 4-1/2" deep x 4" high, making it compact yet spacious enough for essential items.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)