

MPD Tilt Bins Cart Accessory



Stock Photo - Options May Vary

FACTORY NEW

MPD Tilt Bins Cart Accessory. The MPD Tilt Bins Cart Accessory offers a flexible storage solution designed for select MPD medical carts. Each row of tilt bins, displayed on the far left, is sold separately.

These bins provide quick access to frequently used supplies, keeping them visible and organized. Select the appropriate tilt bin style for your specific cart to ensure a proper fit and optimal performance.

STB Tilt Bins fit select carts only and require the separately sold STR-1 Trellis System. Furthermore, for standard carts, use TTB Tilt Bins along with the TTR-6 Trellis System.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)