

MPD Malignant Hyperthermia Cart | MMH-230-DB



Stock Photo - Options May Vary

FACTORY NEW

MPD Malignant Hyperthermia Cart: MMH-230-DB. The MPD MMH-230-DB enhances patient safety by offering a reliable solution to store and transport refrigerated cold packs and saline solutions. This cart keeps emergency supplies at the right temperature and easily accessible during malignant hyperthermia crises.

Its design supports fast response times, helping medical teams act quickly when every second counts. The MMH-230 cart provides ample storage for essential items required in malignant hyperthermia emergencies.

It is built to keep medications and supplies organized and protected, reducing delays during critical care. Durable and easy to move, this cart fits hospital, surgical, and emergency settings for efficient workflow.

FEATURES

- Drawer Configuration: 1-3",1-6"
- 1.6 cubic feet Refrigerator
- Temperature Indicator
- Hospital Grade Plug
- Plastic Top

SPECIFICATIONS

- Drawer Configuration: 1-3",1-6"
- Refrigerator Capacity: 1.6 cubic feet
- Overall Dimensions: 45"H x 24 1/2"D x 29"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)