

MPD Malignant Hyperthermia Cart | MMH-330-DB



Stock Photo - Options May Vary

FACTORY NEW

MPD Malignant Hyperthermia Cart: MMH-330-DB. The MPD MMH-330-DB improves patient safety by keeping cold packs and saline solutions stored and ready for emergencies. This cart allows quick access to critical items during a malignant hyperthermia crisis.

Medical teams can respond faster because refrigerated supplies remain at the correct temperature. Its reliable design ensures that essential items stay organized and protected in high-pressure situations.

The MMH-330 features a 1.6 cubic foot refrigerator/freezer (120V, 60Hz, latex-free) with a temperature indicator for accurate monitoring. It comes with cold packs, a hospital-grade plug with a 10-foot cord, and a security sealing system.

FEATURES

- 1.6 cubic feet Refrigerator/Freezer (120V, 60Hz, Latex-Free)
- Cold Packs
- Temperature Indicator
- Hospital Grade Plug w/10' Cord
- Security Sealing System

SPECIFICATIONS

- 1.6 cubic feet Refrigerator/Freezer: 1.6 cubic feet
- Required Voltage: 120V, 60Hz
- Latex Free
- Drawer Configuration: 2-3" & 1-6"
- Inside drawer dimensions: 22-7/8"D x 16-1/4"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)