

MPD Mini Tower 3 Drawer Isolation Cart | MAT-324



Stock Photo - Options May Vary

FACTORY NEW

The MPD MAT-324 Mini Tower 3 Drawer Isolation Cart provides substantial storage capacity in a compact footprint. This cart features a three-drawer configuration with one 6-inch drawer and two 9-inch drawers, offering deep storage for medical supplies and equipment. The unit is constructed with either aluminum or steel and includes a free extra-thick plastic top with double-wall construction for durability.

The MAT-324 is equipped with premium 5-inch full swivel casters, including two locking casters and one directional caster for controlled mobility. Quiet ball bearing drawer slides ensure smooth operation, while an outrigger stabilizing frame provides structural support. The cart is latex-free and made in the USA.

Available in 17 solid or two-tone color options, the MAT-324 comes with a limited 5-year warranty. The compact dimensions make it ideal for isolation rooms and other space-constrained healthcare environments.

FEATURES

- 3 Drawer Configuration
- Deep Drawers for Maximum Storage
- 5" Swivel Casters (2 casters locking and 1 directional)
- Quiet Ball Bearing Drawer Slides
- Made In The USA

SPECIFICATIONS

- Height: 37 3/4"
- Width: 24 1/2"
- Depth: 24 3/4"
- Drawer: 24"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)