

MPD Mini Tower 4 Drawer Isolation Cart | MAT-430



Stock Photo - Options May Vary

FACTORY NEW

The MPD MAT-430 Mini Tower is a compact 4-drawer isolation cart designed for efficient storage in space-constrained medical environments. It features durable double wall construction, quiet ball-bearing drawer slides for smooth operation, and fully extending drawers for easy access to supplies. The cart includes a free extra-thick plastic top and is available in aluminum or steel construction.

Mobility and safety are enhanced by large 5" premium casters with two locking and one directional caster, plus an outrigger stabilizing frame. The cart offers 30" of total drawer space with a 1-3" and 3-9" drawer configuration, providing flexible storage options for medical supplies and equipment.

The MAT-430 is made in the USA and backed by a five-year limited warranty. Drawer configurations can be customized at no additional charge, and the cart is latex-free for use in sensitive medical settings.

FEATURES

- 4 Drawer Configuration
- Outrigger Stabilizing Frame
- 5" Swivel Casters (2 casters locking and 1 directional)
- Quiet Ball Bearing Drawer Slides
- Made in the USA

SPECIFICATIONS

- Drawer Space: 30"
- Height: 43 3/4"
- Depth: 24 3/4"
- Width: 24 1/2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)