

# MPD Standard 5 Drawer Emergency Cart w/9" Panel | MBS-527P



Stock Photo - Options May Vary

## FACTORY NEW

MPD MBS-527P Standard 5 Drawer Emergency Cart w/9" Panel. MPD Standard Series Carts are full feature models that Healthcare Professionals can use with confidence for years to come. Additionally, All models are available in E-Z Push Aluminum or Steel Construction.

You can customize the drawer configuration at no extra charge. Furthermore, MPD proudly makes all carts in the USA. The MPD MBS-527P features durable steel construction and five locking options.

It includes double wall construction, full-extension ball bearing slides, 5" premium swivel casters (two locking, one directional), and a 5-drawer layout with a 9" panel. Backed by a 5-year limited warranty, it's latex-free and finished in premium powder coat paint.

## FEATURES

- Custom built-to-order
- Steel or Aluminum Construction
- Custom Drawer Configuration
- 9" Panel
- Double Wall Construction

## SPECIFICATIONS

- Drawer Configuration: Four-3", One-6", and One-9"Panel
- Drawer Space: 27"
- Dimensions: 42 3/4"H x 25"D x 32"W

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)