

MPD Standard 6 Drawer Emergency Cart | MBS-630

FACTORY NEW

MPD MBS-630 Standard 6 Drawer Emergency Cart. MPD Standard Series Carts are full feature models that Healthcare Professionals can use with confidence for years to come.

Additionally, all models are available in E-Z Push Aluminum or Steel Construction. Additionally, custom Drawer Configuration can be done without charge and is Proudly Made in the USA.

The MPD MBS-630 Features: Custom Built to Order, Five Locking Options Available, 5-year limited warranty, Double Wall Construction, Full Extension Ball Bearing Drawer Slides, Premium 5" Swivel Casters (2 casters locking and 1 directional), Premium Powder Coat Paint, and is Latex Free. Additionally, the cart is available in your choice of 17 Colors, Solid or Two-Tone.



Stock Photo - Options May Vary

FEATURES

- Custom built-to-order
- 5 Locking Options
- Steel or Aluminum Construction Available
- Double Wall Construction
- Full Extension Ball Bearing Drawer Slides

SPECIFICATIONS

- Dimensions: 39-3/4"H x 25"D x 32"W
- Drawer Configuration: Four-3", One-6" and One-12"
- Drawer Space: 30"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)