

# MPD Standard 4 Drawer Emergency Crash Cart | MBT-424

## FACTORY NEW

The MPD MBT-424 is a 4-drawer emergency crash cart built with either aluminum or steel construction for durability and reliability in healthcare settings. The cart features double wall construction, full extension ball bearing drawer slides, and premium 5-inch swivel casters with two locking and one directional caster for controlled mobility.

The MBT-424 offers five locking options including key lock, breakaway, push-button, electronic, and proximity reader to meet various security requirements. Custom drawer configurations are available to suit specific departmental needs, and the cart is finished with premium powder coat paint in 17 color options, available in solid or two-tone combinations.

This cart is latex-free and backed by a 5-year limited warranty, providing healthcare professionals with a dependable equipment storage solution for emergency and critical care environments.



Stock Photo - Options May Vary

## FEATURES

- 5 Locking Options: Key Lock, Breakaway, Push-Button, Electronic, or Prox Reader
- Double Wall Construction
- Full Extension Ball Bearing Drawer Slides
- Premium 5" Swivel Casters (2 casters locking, 1 caster directional)
- Custom Drawer Configuration available

## SPECIFICATIONS

- Dimensions: 39-3/4"H x 25"D x 32"W
- Drawer Configuration: Two-3", One-6", and One-12"
- Drawer Space: 24"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)