

MPD Standard 5 Drawer Emergency Cart | MBT-524



Stock Photo - Options May Vary

FACTORY NEW

MPD MBT-524 Standard 5 Drawer Emergency Cart. MPD Standard Series Carts are full feature models that Healthcare Professionals can use with confidence for years to come. Additionally, all models are available in E-Z Push Aluminum or Steel Construction.

You can configure your custom drawer on the MBT-524 at no extra charge. In addition, the cart is available in 17 Colors, Solid or Two-Tone. Furthermore, see below for additional color information.

The MPD MBT-524 Features: 5 Locking Options Available, comes with a 5-year limited warranty, Double Wall Construction, Full Extension Ball Bearing Drawer Slides, Premium 5" Swivel Casters (2 casters locking and 1 directional), Premium Powder Coat Paint, and Latex Free.

FEATURES

- Custom built-to-order
- Five Locking Options
- Steel or Aluminum Construction Available
- Double Wall Construction
- Full Extension Ball Bearing Drawer Slides

SPECIFICATIONS

- Drawer Configuration: Three-3", One-6", and One-9"
- Drawer Space: 24"
- Dimensions: 39-3/4"H x 25"D x 32"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)