

# MPD Mini Tower 5 Drawer Emergency Cart w/Breakaway Lock | MCT-524



Stock Photo - Options May Vary

## FACTORY NEW

MPD MCT-524 5 Drawer Emergency Cart w/ Lock. MPD has manufactured High-Quality Medical Carts and Accessories for healthcare professionals for over twenty Years.

MPD's Mini Carts offer the best solution to space limitations with a slimmer profile and compact 5" Full Swivel Casters (two casters are locking and one-directional). The MPD MCT-524 Features Double-Wall Construction, a Small footprint, five Drawers, an Outrigger Stabilizing Frame, and a Breakaway Lock.

All Mini Carts are available in E-Z Push Aluminum or Steel construction and a range of drawer configurations. In addition, custom drawer configurations can be done without charge.

## FEATURES

- Custom built-to-order
- Choose from 17 Colors
- Steel or Aluminum Construction Available
- Double Wall Construction
- Full Extension Ball Bearing Drawer Slides

## SPECIFICATIONS

- Drawer Configuration: Three-3", One-6", and One-9"
- Drawer Space: 24"
- Dimensions: 37-3/4"H x 24-3/4"D x 24-1/2"W

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)