

MPD Mini Tower 6 Drawer Cart with Break-Away Lock | MCT-630



Stock Photo - Options May Vary

FACTORY NEW

MPD 6 Drawer Cart w/ Lock: MCT-630. MPD Mini Carts provide the best solution for space limitations with a slim profile and 3" full swivel casters. Mini Carts are available in E-Z Push Aluminum or Steel construction with a range of drawer configurations.

Choice of 17 Colors, Solid or Two-Tone. The MPD MCT-630 Features a breakaway Lock, Quiet Ball ball-bearing drawer Slides offering a Full Extension Drawer Opening, a Stabilizer Frame, Full Swivel Casters, and a FREE Extra Thick Plastic Top.

Furthermore, Custom Drawer configuration can be done without charge. Furthermore, we proudly make Mini Carts in the USA.

FEATURES

- Proudly Made in the USA
- Steel and E-Z Push Aluminum Construction Available
- Quiet Ball Bearing Slides
- Full Drawer Extensions
- Stabilizer Frame

SPECIFICATIONS

- Drawer Space: 30"
- Height: 43-3/4"
- Depth: 24-3/4"
- Width: 24-1/2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)