

# MPD Mini 5 Drawer Emergency Cart w/Breakaway Lock | MCC-524



Stock Photo - Options May Vary

## FACTORY NEW

MPD MCC-524 Mini 5-Drawer Cart w/ Lock. MPD has manufactured High-Quality Medical Carts and Accessories for healthcare professionals for over twenty Years.

MPD's Mini Carts offer the best solution to space limitations with a slimmer profile and compact 3" Full Swivel Casters (2 casters locking). Additionally, All MPD Carts are proudly made in the USA.

The MPD MCC-524 Features Double Wall Construction, a Small Footprint, 5 Drawers, and a Breakaway Lock. Mini Carts come in E-Z Push Aluminum or Steel construction and offer a variety of drawer configurations. Furthermore, you can request a custom drawer configuration at no additional cost.

## FEATURES

- Custom built-to-order
- Choose from 17 Colors
- Breakaway Lock
- Steel or Aluminum Construction Available
- Double Wall Construction

## SPECIFICATIONS

- Drawer Configuration: Three-3", One-6", One-9"
- Drawer Space: 24"
- Dimensions: 34"H x 18"D x 18-1/2"W
- 5-Year Limited Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)