

MPD Standard 4 Drawer Isolation Cart | MKS-430-Y



Stock Photo - Options May Vary

FACTORY NEW

MPD MKS-430-Y Standard 4 Drawer Isolation Cart. MPD Standard 4 Drawer Isolation Cart: Standard Series Carts are full-feature models that Healthcare professionals can use with confidence for years to come. All models are available in E-Z Push Aluminum or Steel construction.

Additionally, the cart has an extra thick plastic top and Double-Wall Construction, with Premium 5" full-swivel casters; two are casters locking, and one is directional. Additional MPD MKS-430-Y features an Outrigger Stabilizing Frame and Quiet Ball-bearing Drawer Slides.

Furthermore, the cart displays Optional Accessories: a Coat Rack, Railing, and Utility Hooks, sold separately and not included. Additionally, It is proudly made in the U.S.A.

FEATURES

- 4 Drawer Configuration
- The cart is Latex Free
- Deep Drawers for Maximum Storage
- E-Z Push Aluminum or Steel Construction
- Steel Double Wall Construction

SPECIFICATIONS

- Drawer Configuration: One-3" and Three-9"
- Drawer Space: 30"
- Dimensions: 46-1/2"H x 25"D x 32"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)