

# MPD Standard 4 Drawer Medication Cart | MKS-430TL



Stock Photo - Options May Vary

## FACTORY NEW

MPD MKS-430TL Standard 4 Drawer Medication Cart. Standard Series Carts are Full Feature Models that Healthcare Professionals can use with confidence for years to come.

Furthermore, All models are available in EZ-Push Aluminum (price reflects Aluminum Cart) or Steel Construction (see options). Additionally, All MPD Standard Carts are Proudly Made in the USA and are totally latex-free. The MPD MKS-430TL Features: The cart features a free extra-thick plastic top, custom drawer configuration, and double-wall construction.

It includes twin poly handles, premium 5" full-swivel casters (2 with toe-locking brakes), an outrigger stabilizing frame, and quiet ball-bearing drawer slides with full extension.

## FEATURES

- Made In The USA
- Latex Free
- Aluminum Construction
- Quiet Ball Bearing Slides
- Full Drawer Extensions

## SPECIFICATIONS

- Drawer Space: 30"
- Height: 46"
- Depth: 25-1/2"
- Width: 32"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)