

MPD Mini Isolation 3 Drawer Cart | MMC-324



Stock Photo - Options May Vary

FACTORY NEW

MPD MMC-324 Mini Isolation 3 Drawer Cart. The MPD MMC-324 offers exceptional functionality for healthcare professionals. Part of the MPD Standard Series, this cart provides years of reliable use.

It comes in both aluminum and steel construction for durability, and you can choose from 17 solid or two-tone colors. Key features of the MMC-324 include a free, extra-thick plastic top and double-wall construction for added strength.

It also has premium 5" full swivel casters, two of which lock, and one that is directional for easy maneuvering. The outrigger stabilizing frame ensures the cart remains steady, while the quiet ball-bearing drawer slides allow for smooth operation.

FEATURES

- 3 Drawers
- Drawer Configuration: 1-6" and 2-9"
- The cart is Latex Free
- Deep Drawers for Maximum Storage
- E-Z Push Aluminum or Steel

SPECIFICATIONS

- Drawer Space: 24"
- Height: 34"
- Depth: 18"
- Width: 18 1/2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)