

MPD Standard Unit Dose Medication Cart Package | MUS-330-BPKG



Stock Photo - Options May Vary

FACTORY NEW

MPD MUS-330-BPKG Standard Unit Dose Medication Cart Package. Unit Dose Meets Economy. Standard Series Carts are Full-Feature Models that Healthcare Professionals can use with confidence for years to come.

All models are available in EZ-Push Aluminum or Steel Construction. Moreover, all MPD carts are proudly made in the USA and are completely latex-free. The MPD MUS-330-BPKG Features: Carts feature steel and E-Z Push aluminum construction with ball bearing slides and full drawer extensions.

They include stabilizer frames, full-swivel casters, and extra-thick plastic tops. A pull-out cassette shelf allows easy access and exchange. Additionally, the package includes a cart, cassette, 20 bins, and a 3" lockable drawer.

FEATURES

- Made In The USA
- Latex Free
- Includes Cart
- Steel and Aluminum Construction
- Quiet Ball Bearing Slides

SPECIFICATIONS

- Height: 46-1/4"
- Depth: 25"
- Width: 32"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)