

MPD Value 4 Drawer Cart | MVP-424-B



Stock Photo - Options May Vary

FACTORY NEW

MPD MVP-424-B Value 4 Drawer Cart. MPD has over 20 Years of Experience Manufacturing Premium Quality Medical Carts and Workstations including Super Saver Carts for budget-conscious facilities, along with a wide variety of Cart Accessories to maximize storage and organization.

Currently, there are over 2000,00 MPD products in use today in hospitals and medical facilities coast-to-coast and around the world. The MPD MVP-424-B Features: Strength and Versatility in every cart, Steel Construction, Includes Plastic Top, Locking Top Drawer, Four Drawers (lock secures Top Drawer Only, all other drawers are Non-Locking), and 24" of Drawer Space.

Additionally, available in Beige Only. Note: The MVP-424-B is shown with a Microscope, sold separately.

FEATURES

- Proudly Made in the USA
- Steel Construction
- Includes Plastic Top
- Locking Top Drawer
- Drawer Configuration: One-3", Two-6", and One-9"

SPECIFICATIONS

- Height: 34-1/2"
- Depth: 18"
- Width: 18-1/2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)