

MPD Select 5 Drawer Aluminum Emergency Cart | SBT-524A



Stock Photo - Options May Vary

FACTORY NEW

MPD Select 5 Drawer Emergency Cart: SBT-524A. MPD has manufactured High-Quality Medical Carts and Accessories for healthcare professionals for over 20 Years. MPD proudly makes all carts in the USA. The following emergency/code/crash cart is no exception! Strength and Versatility in every cart.

Additionally, a Five-Year Limited Warranty is included. The MPD SBT-524A Features: 5 Drawers to Organize Supplies and six Locking Options Available (Key Lock, Push-Button, Breakaway, Key Lock Bars, Electronic, or Prox Reader with Key Override). Twin Poly Push Handles, Plastic Top, Dual Pull Out Side Shelf, and Full Extension Ball Bearing Drawer Slides.

Furthermore, all Select Carts are available in 17 Solid or Two-Tone Colors.

FEATURES

- Proudly made in the USA
- Custom Built to Order
- Available in 17 Colors
- Six Locking Options
- Steel or Aluminum Construction

SPECIFICATIONS

- Drawer Configuration 3-3" 1-6" 1-9"
- Dimensions: 40-5/8"H x 25"D x 32"W
- 5 Year Limited Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)