

MPD Select 6 Drawer Aluminum Emergency Cart | SBT-627A



Stock Photo - Options May Vary

FACTORY NEW

MPD Select 6 Drawer Emergency Cart: SBT-627A. MPD has manufactured High-Quality Medical Carts and Accessories for healthcare professionals for over 20 Years. All carts are Proudly Made in the USA. The following emergency/code/crash cart is no exception.

The MPD SBT-627A Features: Six Locking Options Available (Key Lock, Push-Button, Breakaway, Key Lock Bars, Electronic, and Prox Reader with Key Override). Twin Poly Push Handles, Plastic Top, Dual Pull Out Side Shelf, and Full Extension Ball Bearing Drawer Slides.

FEATURES

- All MPD Carts are Proudly Made in the USA
- Six Locking Options Available
- Steel or aluminum Construction Available
- Double Wall Construction
- Full Extension Ball Bearing Drawer Slides

SPECIFICATIONS

- Drawer Configuration: Four-3", One-6", and One-9"
- Dimensions: 44"H x 25"D x 32"W
- 5 Year Limited Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)