

# MPD Select 5 Drawer Aluminum Cart w/Electronic Lock | SET-524A



Stock Photo - Options May Vary

## FACTORY NEW

MPD Select 5 Drawer Aluminum Cart w/Electronic Lock: SET-524A. With 20+ years of experience, MPD delivers high-quality medical carts and accessories.

The premium Select Series meets top industry standards. Moreover, proudly made in the USA, this cart features expert craftsmanship, dual pull-out side shelves, and ball-bearing slides for smooth, functional use in clinical settings.

SET-524A Features: This cart includes 5 drawers with 24" of storage space to organize a wide range of supplies. Furthermore, the electronic lock offers dual-code access, battery operation with no recharging, and a low battery indicator for convenience.

## FEATURES

- Made in the USA
- Latex-Free
- Electronic Lock
- Multi-code Capacity
- Battery Operated, No Recharging Needed

## SPECIFICATIONS

- Drawer Configuration: Three-3", One-6", and One-9"
- Drawer Space: 24"
- Cart Dimensions: 40-5/8"H x 25"D x 32"W

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)