

MPD Select 5 Drawer Aluminum Mobile Anesthesia Workstation | SKT-524A



Stock Photo - Options May Vary

FACTORY NEW

MPD Select 5 Drawer Aluminum Mobile Anesthesia Workstation: SKT-524A. MPD has over 20 Years of Experience manufacturing premium quality medical carts, narcotic cabinets, mobile workstations, and accessories. The Select Series is MPD's newest Premium Line of Carts and Accessories.

SKT-524A Features Key Lock, 5 Drawers with 24" Drawer Space. Also features Deluxe Drawer Pulls, each with a convenient area for a label and a large area at the top front for identifying the cart.

Select Carts have Dual Pull Out Side Shelves with Quiet Ball Bearing Slides. Furthermore, features Three Adjustable Side Accessory Panels for Attaching Accessories. Customize your Drawer Configurations to meet your individual needs is available.

FEATURES

- Made in the USA
- Latex-Free
- User-Friendly
- E-Z Push Aluminum Construction
- 5 Drawers

SPECIFICATIONS

- Drawer Configuration: Three-3", One-6", and One-9"
- Drawer Space: 24"
- Cart Dimensions: 40 5/8"H x 25"D x 32"W
- 5 Drawer Key Lock

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)