

# MPD Select 6 Drawer AL Anesthesia Cart w/Electronic Lock | SES-627A



Stock Photo - Options May Vary

## FACTORY NEW

MPD Select 6 Drawer AL Anesthesia Cart w/Electronic Lock: SES-627A. MPD's Select Series represents over 20 years of expertise in crafting premium medical carts and accessories. The SES-627A anesthesia cart offers top-tier performance with a sleek design.

It features dual pull-out side shelves, glide smoothly on quiet ball-bearing slides. Healthcare providers rely on MPD for durable, functional solutions that improve workflow and meet demanding clinical environments. This anesthesia cart includes an advanced electronic lock with key override, multi-code access, and a tamper indicator.

It runs on batteries and needs no recharging. Additionally, it auto-locks and shows a low battery warning.

## FEATURES

- Made in the USA
- Latex-Free
- Electronic Lock
- Multi-code Capacity
- Battery Operated, No Recharging Needed

## SPECIFICATIONS

- Drawer Configuration: Four - 3", One - 6", and One - 9"
- 27" Drawer Space
- 44"H x 25"D x 32"W

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)