

# MPD Select 6 Drawer AL Anesthesia Cart w/Electronic Lock | SET-630A



Stock Photo - Options May Vary

## FACTORY NEW

MPD Select 6 Drawer AL Anesthesia Cart w/Electronic Lock: SET-630A. MPD's Select Series combines over 20 years of medical cart expertise with advanced design and functionality. The SET-630A anesthesia cart includes dual pull-out side shelves on quiet ball-bearing slides.

Moreover, proudly made in the USA, this premium cart improves efficiency in fast-paced medical environments. It offers a perfect blend of durability, organization, and smooth mobility for healthcare professionals.

The SET-630A comes with a secure electronic lock and a free plastic top for extra workspace. Additionally, it has six roomy drawers with 30 inches of storage space.

## FEATURES

- Made in the USA
- Latex-Free
- Electronic Lock
- Free Plastic Top
- Dual Pull Out Side Shelves

## SPECIFICATIONS

- Drawer Configuration: Three 3", Two -6", and One 9"
- 30" Drawer Space
- Cart Dimensions: 47 1/4"H x 25"D x 32"W

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)