

MPD Select 6 Drawer Aluminum Anesthesia Cart w/ Push Button Lock | SIS-627A



Stock Photo - Options May Vary

FACTORY NEW

MPD Select 6 Drawer Aluminum Anesthesia Cart w/ Push Button Lock: SIS-627A. MPD has over 20 Years of Experience manufacturing premium quality medical carts, narcotic cabinets, mobile workstations, and accessories. The Select Series is MPD's newest Premium Line of Carts and Accessories.

This cart is proudly made in the USA and is Latex-Free. MPD's Select Carts offer the very best and most user-friendly carts. Features Include: User-Friendly Push-Button Locking - Easy and convenient mechanical push-button locking with key lock override.

Simply enter the combination, turn the knob one-quarter turn and your cart is open, one quarter-turn back and your cart is locked. Best of all...NO Batteries!

FEATURES

- Made in the USA
- Latex-Free
- User-Friendly
- E-Z Push Aluminum Construction
- 6 Drawers

SPECIFICATIONS

- Drawer Configuration: Four-3", One-6", and One-9"
- Drawer Space: 27"
- Cart Dimensions: 44 1/4"H x 25"D x 32"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)