

MPD Standard 4 Drawer Anesthesia Cart w/ Electronic Lock | MET-430



Stock Photo - Options May Vary

FACTORY NEW

The MPD MET-430 is a 4-drawer anesthesia cart with durable steel and aluminum construction designed for secure medication storage in healthcare settings. It features an electronic lock with multi-code capacity, dual-code programming, and a tamper indicator. The battery-operated lock requires no recharging and includes an auto-lock system with backup key access.

The cart offers full-extension drawers for quick access to supplies and is built with a stabilizer frame and full-swivel casters for smooth movement and stability. The extra-thick plastic top withstands daily use in fast-paced environments. Quiet ball bearing slides ensure smooth drawer operation.

Manufactured in the USA with 100% latex-free materials, the MET-430 is suitable for allergy-sensitive areas. It is available in 17 solid or two-tone colors to match various medical environments.

FEATURES

- Electronic Lock
- Made in the USA
- Latex Free
- Full Drawer Extensions
- Full-Swivel Casters

SPECIFICATIONS

- Drawer Configuration: One-3" and Three-9"
- Drawer Space: 30"
- Cart Dimensions: 46-1/4"H x 25"D x 32"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)