

MPD Standard 5 Drawer Anesthesia Cart w/ Push Button Lock | MIS-524



Stock Photo - Options May Vary

FACTORY NEW

MPD Standard 5 Drawer Anesthesia Cart w/ Push Button Lock: MIS-524. The MPD Standard MIS-524 delivers strength, security, and long-term performance. It features full-extension drawers on quiet ball-bearing slides, allowing quick access to supplies. The mechanical push-button lock with key override ensures reliable security without batteries.

Enter the code, turn the knob, and you're ready to go. This cart suits fast-paced healthcare environments that demand dependability and ease of use. The MIS-524 is designed for easy movement and efficiency.

It has five drawers, a stabilizer frame, and wheels that turn smoothly for easy transport. You can choose a strong steel frame or a lighter EZ-Push aluminum option.

FEATURES

- 5 Drawers
- Made in the USA
- Latex Free
- Steel and Aluminum Construction
- Quiet Ball Bearing Slides

SPECIFICATIONS

- Drawer Configuration: Three-3", One-6", and One-9"
- Drawer Space: 24"
- Cart Dimensions: 39-3/4"H x 25"D x 32"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)