

MPD Standard 6 Drawer Anesthesia Cart w/ Key Lock | MKT-627



Stock Photo - Options May Vary

FACTORY NEW

MPD Standard 6 Drawer Anesthesia Cart w/ Key Lock: MKT-627. The MPD MKT-627 delivers durability and efficiency for hospitals, clinics, and long-term care settings. This full-feature cart is available in EZ-Push aluminum or steel construction to fit different clinical needs. With various heights and customizable drawer configurations, the MKT-627 supports reliable daily use and secure medical storage.

Key features include six smooth-gliding, full-extension drawers and a customizable layout. Stabilizer frames and swivel casters ensure smooth transport. Furthermore, the cart includes a durable, easy-to-clean plastic top.

Additionally, a secure key lock system makes it ideal for medication management and anesthesia storage.

FEATURES

- Made in the USA
- Latex Free
- Steel and Aluminum Construction
- Quiet Ball Bearing Slides
- Full Drawer Extensions

SPECIFICATIONS

- Drawer Configuration: Four-3", One-6", and One-9"
- Drawer Space: 27"
- Cart Dimensions: 42 3/4"H x 25"D x 32"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)