

# MPD Standard 6 Drawer Anesthesia Cart w/ Proximity Reader | MXT-630



Stock Photo - Options May Vary

## FACTORY NEW

MPD Standard 6 Drawer Anesthesia Cart w/ Proximity Reader: MXT-630. Standard Series Carts are Full Feature Models that Healthcare Professionals can use with confidence for years to come.

All models are available in EZ-Push Aluminum or Steel Construction with a wide range of heights, and drawer configurations. The MXT-630 Features: Electronic Lock, EZ Open Auto Lock w/Key Override, Free Extra Thick Plastic Top, Custom Drawer Configuration can be done without charge and cart has double-wall construction, Twin Poly Handles, and Premium 5" Full Swivel Casters (2 locking casters, 1 directional).

Outrigger Stabilizing Frame and Quiet Ball Bearing Drawer Slides offer a Full Extension Drawer Opening. Note: Proximity Reader Software and Card Sold Separately.

## FEATURES

- Made in the USA
- Latex Free
- Proximity Cart
- 6 Drawers

## SPECIFICATIONS

- Drawer Configuration: Three-3", Two-6", and One-9"
- Drawer Space: 30"
- Dimensions: 46-1/4"H x 25"D x 32"W

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)