

MPD Standard 6 Drawer Anesthesia Cart w/ Push Button Lock | MIT-627



Stock Photo - Options May Vary

FACTORY NEW

The MPD Standard 6 Drawer Anesthesia Cart is a full-featured mobile workstation designed for healthcare professionals. Built with choice of steel or aluminum construction, the cart features six drawers with quiet ball bearing slides and full extensions for convenient access to supplies and equipment.

The cart includes a mechanical push-button locking system with key lock override, requiring no batteries for operation. A simple combination entry followed by a quarter-turn of the knob opens the cart, with another quarter-turn to secure it.

Standard features include a stabilizer frame, full-swivel casters, and an extra-thick plastic top for durability. The cart is made in the USA, latex-free, and available in 17 solid or two-tone color options.

FEATURES

- Made in the USA
- Latex Free
- Quiet Ball Bearing Slides
- Full Drawer Extensions
- Push Button Lock

SPECIFICATIONS

- Drawer Configuration: Four-3", One-6", and One-9"
- Drawer Space: 27"
- Cart Dimensions: 42-3/4"H x 25"D x 32"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)