

MPD Standard 6 Drawer Anesthesia Cart w/ Push Button Lock | MIT-630



Stock Photo - Options May Vary

FACTORY NEW

The MPD Standard 6 Drawer Anesthesia Cart (MIT-630) is designed for reliable performance in demanding clinical environments. This cart features a push-button lock for secure access, double-wall construction, and an extra-thick plastic top for durability. The unit is available in either steel or aluminum construction with a choice of 17 solid or two-tone colors.

The MIT-630 includes six drawers with quiet ball bearing slides that extend fully for easy access to supplies. Each drawer can be customized without additional charge to meet your specific organizational needs. The cart is equipped with premium 5-inch full swivel casters, including two locking casters and one directional caster for optimal mobility and stability.

Additional features include twin poly handles, an outrigger stabilizing frame, and a compact footprint suitable for various clinical settings. The cart's robust construction and thoughtful design make it a dependable solution for anesthesia supply management and organization.

FEATURES

- Push Button Lock
- Quiet Ball Bearing Slides
- Full Drawer Extensions
- Stabilizer Frame
- Extra Thick Plastic Top

SPECIFICATIONS

- Dimensions: 46-1/4"H x 25"D x 32"W
- Drawer Configuration: Three-3", Two-6", and One-9"
- Drawer Space: 30"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)