

# MPD Standard Wide Aluminum Punch-Card Medication Cart | WKT-430A



Stock Photo - Options May Vary

## FACTORY NEW

MPD WKT-430A Standard Wide Aluminum Punch-Card Medication Cart. Maximize Your Storage with a Super Wide MPD Cart. All MPD Carts are Proudly Made in the USA and are totally latex-free.

Carts feature steel and E-Z Push aluminum construction with quiet ball-bearing slides and full drawer extensions. Custom Drawer Configuration is Available, The MPD WKT-430A Features Standard Wide carts with stabilizer frames and full-swivel casters, as well as FREE extra-thick plastic tops.

They also have Twin Poly Push Handles, an Outrigger Stabilizing Frame, and Double-Wall Construction. Additional Features: This Wide Aluminum Punch-Card Medication Cart fits 12 optional TPC-1 punch-card racks that fit up to 25-30 punch cards per rack.

## FEATURES

- Made In The USA
- Latex Free
- Steel and Aluminum Construction
- Quiet Ball Bearing Slides
- Full Drawer Extensions

## SPECIFICATIONS

- Drawer Space: 30"
- Height: 46-1/4"
- Depth: 25-1/2"
- Width: 40-1/2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)