

Oakworks Lithotripsy-Urology Imaging Table: CFLU401

FACTORY NEW

The Oakworks CFLU401 is a versatile lithotripsy-urology imaging table designed for both mobile and stationary use. It features four powered motions including height adjustment, Trendelenburg tilt, longitudinal travel, and lateral travel, with a 500 lbs patient weight capacity and height range of 27.5" to 45.5". The table's easy-to-use extension system accommodates oversized or undersized patients and enables left or right-sided procedures.



The CFLU401 quickly converts to a fully-featured urology table for cystoscopies and other urology procedures. It includes a large carbon fiber top, hand rail, lightweight end extensions, six radiolucent extensions (each supporting 200 lbs), T-rail slide, lock lateral extensions, and a digital angle gauge for precise positioning.

The table offers convenient hand and foot controls for all functions and includes a 1" Comfort Foam table top pad, urology bag hoop, patient safety strap, and heavy-duty 5" casters. It is FDA listed, cETL approved, and CE marked.

Stock Photo - Options May Vary

FEATURES

- 4 Powered Motions
- Motions Included: Height, Trendelenburg tilt, Longitudinal travel, Lateral travel
- Patient Weight Capacity: 500 lbs.
- 27.5" to 45.5" Height Range
- All Motion Hand and Foot Control
- Large Carbon Fiber Top
- Digital Angle Gauge accurately displays the angle
- Easily Converts to a fully-featured Urology table in seconds

SPECIFICATIONS

- Carbon Fiber Top: 23" wide x 47" Long
- Imaging Area: 19" Wide x 23.5" Long
- Patient Weight Capacity: 500 lbs - Table maximum, with or without extensions
- Each extension: 200 lbs.
- Table Weight: 618 lbs
- Electrical: 120V, 8.2 amps, 60Hz operation
- Warranty: 3 Years Parts, 2 Year Labor
- Certifications: FDA Listed, cETL Approved, CE Marked

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com