

# Oakworks Aesthetic Procedure Chair | 100A Series

## FACTORY NEW

The Oakworks 100A Series Procedure Chair combines medical functionality with patient comfort through spa-like design and ergonomic engineering. The chair features powered height adjustment, Fowler positioning, and Trendelenburg tilt capabilities, all controlled via multi-function foot pedal for hands-free operation. An open base design supports proper provider ergonomics while the acrylic base cover and skirt allow customization to match facility décor.



Stock Photo - Options May Vary

The 100A includes a full 74-inch exam table with adjustable headrest, manual leg rest, and armrests for patient support during supine or prone procedures. The chair accommodates patients up to 600 pounds and features 2.5-inch thick multi-layered padding upholstered in premium, PVC-free TerraTouch material available in 21 colors.

Powered adjustments include an 80-degree backrest range and 15-degree Trendelenburg with 5-degree reverse Trendelenburg capability. The chair operates on standard 115V/60Hz power and meets regulatory certifications including ETL testing to UL 60601-1 and CE marking.

## FEATURES

- 100A Series Power Procedure Chair with \*Acrylic Upgrade Included
- Powered Height, Fowler & Trendelenburg
- Multi-function Foot Control
- Full 74" Long Exam Table
- Available in 21 Tera Touch™ Upholstery Colors

## SPECIFICATIONS

- 19" - 35" Powered Height Range
- Patient Capacity: 600 lbs
- 3 Section Top: 74"L (not including extended Head Rest) x 27"W
- Padding: 2.5" Thick, Multi-Layered
- Table Voltage: 115V/60Hz

*\*Brochure and images are a general overview of the product. Actual product options may vary.*