

# Oakworks Echocardiography Ultrasound Table

## FACTORY NEW

The Oakworks Echocardiography Table is designed for cardiac and vascular ultrasound procedures with a sliding cardio panel that opens up to 15 inches to accommodate both left and right-handed scanning. The electric height adjustment ranges from 22 to 38 inches, promoting safe wheelchair transfers and easy adaptation between vascular and cardiac applications. A dual tower design provides strength and rigidity with a 550-pound weight capacity.



*Stock Photo - Options May Vary*

The table features a two-section top with electric Fowler positioning controlled via hand control, and adjustable Trendelenburg and reverse tilt capabilities of  $\pm 15$  degrees. The 30-inch width and 72-inch length accommodate most patients, with an optional extender available for taller individuals, extending the table to 84 inches.

The table is FDA-listed, ETL-approved, and CE-marked, with a 3-year warranty. Additional features include locking casters for mobility and a paper roll holder for convenience during procedures.

## FEATURES

- Electric Height Range: 22" - 38"
- Patient Weight Capacity: 550 lbs
- Trendelenburg / Reverse Tilt:  $\pm 15^\circ$
- Electric Fowler Hand Control (All The table movements)
- Sliding Cardio Panel - Manual
- Two Section Top
- Locking Casters
- FDA Listed, ETL Approved, CE Marked

## SPECIFICATIONS

- Patient Weight Capacity: 550 lbs
- Height Range: 22" - 38"
- Width: Only 30"
- Length: 72"
- Top 2 Section: Fowler - 34", Seat - 38"
- Trendelenburg/Reverse;  $\pm 15^\circ$
- P/N: 65258

*\*Brochure and images are a general overview of the product. Actual product options may vary.*