

Oakworks Vascular Table with Fowler



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks Vascular Table with Fowler and Feature Package is designed for vascular and ultrasound imaging procedures. Dual towers and electric Fowler positioning enable up to 70° reverse Trendelenburg positioning for effective standing and full weight-bearing patient positioning. The split drop-away footboard facilitates safe non-weight-bearing scans while the table's height range of 19" to 33" reduces technologist strain during lower extremity vein access.

This table accommodates bariatric patients up to 550 lbs and features electric height adjustment, ±18° Trendelenburg tilt, and up to 40° reverse Trendelenburg tilt. The integrated split footboard, locking casters, dual side rails, and patient restraint straps provide stability and safety during procedures.

The Feature Package includes hand control for all movements, 4-inch locking casters, paper roll holder, and 23.5-inch side rails. The table comes with 3 years of parts coverage and 1 year of labor warranty. It is FDA listed, cETL approved, and CE marked.

FEATURES

- Electric Height Range: 19" - 33"
- Reverse Trendelenburg Tilt: 40° + Electric Fowler = 70°
- Patient Weight Capacity: 550 lbs.
- Integrated footboard (split, drop away)
- Locking Casters: 4" diameter double wheel
- Hand Control (All movements)
- Warranty: 3 Years Parts, 1 Year Labor
- FDA Listed, ETL Approved, CE Marked

SPECIFICATIONS

- Patient Weight Capacity: 550 lbs
- Height Range: 19"- 33"
- Width: 26" or 30"
- Length: 72"
- Trendelenburg: ±18°
- Reverse Trendelenburg: up to 70°
- Padding: 3" multi-layer
- Table Weight: 359 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)