

Oakworks CFPM100 Imaging Table w/Integrated Headrest



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks CFPM100 Imaging Table is a single-tower lift system designed for pain management, surgical centers, and multidisciplinary facilities. It supports a 500-pound patient capacity and is compatible with all mobile C-arm imaging systems. The table features an integrated headrest and electric height adjustment.

The CFPM100 accommodates a height range of 25 inches to 43 inches, enabling ergonomic working positions while facilitating safe wheelchair transfers at its lowest setting. The table includes a soft-touch base material that protects against C-arm damage and scratches during routine use.

Standard inclusions with purchase are a 1-inch Radiolucent Comfort Foam tabletop pad, locking casters, a crescent face rest pad, and two patient safety straps. This system provides a cost-effective solution for new institutions, expanding facilities, and startup practices requiring an elevating imaging table.

FEATURES

- Integrated Headrest Top
- The single-tower lift system
- Electric Height Adjustment with a range of 22" to 40"
- Weight capacity: 500 lbs distributed evenly
- Compatible with all mobile c-arm imaging systems
- Soft-touch table base material
- Durable base material Protects and Minimizes c-arm scrapes and damage
- Comes standard with Foot Control

SPECIFICATIONS

- Carbon Fiber Top with Headrest Dimensions: 22" W x 84" L
- Imaging Space: 56" of Metal-Free Space
- Patient Weight Capacity: 500 lbs distributed evenly
- Electronics: 110 Volt Motor
- 4" Locking Casters: Locking / Unlocking casters for mobile or permanent positioning
- Table Weight: 343 lbs
- Frame: Single Tower Steel Base
- Warranty: 3 years parts, 2 years labor

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)