

Oakworks CFPM300 Imaging Table: Rectangular Top

FACTORY NEW

The Oakworks CFPM300 Imaging Table features a rectangular carbon fiber top with three powered motions: height adjustment, Trendelenburg tilt, and lateral tilt. The table offers an industry-leading 18-inch height range (26-43 inches) and a 56-inch metal-free imaging area, supporting a 500-pound patient weight capacity.

The table includes a hand control for all three motions, 4-inch casters for mobility, and a patient safety strap with Velcro closure. A 1-inch comfort foam pad is provided as standard equipment.

Manufactured in the USA, the CFPM300 is FDA listed, cETL approved, and CE marked. It comes with a 3-year parts and 2-year labor warranty.



Stock Photo - Options May Vary

FEATURES

- 3 Powered Motions: Height, Trendelenburg Tilt, and Lateral Tilt
- Industry's Widest Adjustable Height Range of 18" (26 in - 43 in (66-109 cm))
- Patient Weight Capacity: 500 lbs (227 kg).
- 56" of Metal Free Imaging Area
- Patient Safety Strap with Velcro® Closure
- Hand Control Included
- 4" Casters
- Made in the USA

SPECIFICATIONS

- Patient Weight Capacity: 500 lbs
- Rectangular Carbon Fiber Top: 24" (61 cm) wide x 80" (203 cm) long.
- Movement 1 - Height: 26 in - 43 in (66-109 cm) Height Range without pad.
- Movement 2 - Trendelenburg Tilt: +15°/-12°
- Movement 3 - Lateral Tilt: ±15°
- Controls: 1 Hand (3 Motions)
- Table Weight: 475 lbs (216 kg).
- Warranty: 3 Years Parts, 2 Year Labor

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)