

Oakworks Prone Pillow



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks Prone Pillow provides superior comfort and proper patient positioning for prone or supine procedures. The advanced lower profile design reduces excessive cervical extension and flexion while the resilient platform delivers optimal support for patients of any size or shape. Soft contact surfaces with large vent areas ensure exceptional support and easy breathing during treatment.

FEATURES

- Superior patient positioning for any size/shape patient
- Soft contact surfaces offer exceptional support
- The design has large vent areas for easy breathing
- Can be used with any fluoroscopic or treatment table
- Low profile reduces cervical extension/flexion

SPECIFICATIONS

- Weight: 3 LBS
- Measurements - Platform: 19" x 19" x 1/8 " thick
- Crescent Face Pad: 11" diameter
- Adjuster Pad: 19" x 3" x 6"
- Padding: AeroCel™
- Upholstery: TerraTouch™
- Optional Face Rest Covers Available
- P/N: 66180

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)