

Oakworks Spine Positioning System II



Stock Photo - Options May Vary

FACTORY NEW

The Oakworks Spine Positioning System II is an integral component of pain management suites designed to enhance procedural efficiency and patient outcomes. The system reduces setup time, improves patient comfort, and minimizes unwanted movement during spine procedures. Clearer visualization of target anatomy enables physicians to perform spinal injections more efficiently and securely.

The system enhances visualization of key anatomic landmarks for safe, accurate spinal injections while addressing imaging, stability, and patient comfort issues. The design improves visibility of the cervical, thoracic, and lumbar spine, including the cervicothoracic junction, and ensures optimal positioning regardless of patient body habitus.

The complete 8-piece kit includes a portable case for clean transport and storage. The system is FDA Listed and CE Marked, with a 2-year warranty covering frame, fabric, and padding.

FEATURES

- An integral component of the pain management suite
- Improves visualization of the necessary anatomic landmarks for safe and accurate spinal injections
- Solves imaging, stability, and patient comfort issues associated with spine procedures
- The design provides greater visualization of the cervical, thoracic, and lumbar spine including the cervicothoracic junction regardless of the patient's body habitus
- Radiolucent Frame has one cam lock which facilitates cervical flexion and extension
- Adequate free space under see-through face section for aeration, supplemental O2, and jaw/face
- Prevents patient overheating and enables easier communication with the patient without sound muffling
- 8-piece kit that comes with a portable case for clean transport and storage

SPECIFICATIONS

- Weight: 16 lbs
- Frame with face rest (1): 12" x 32.5"
- Crescent Face Pad (1): 12" diameter
- Contoured Torso Support Pad: 6.5" x 23" x 30"
- Frame: Radiolucent
- Padding: Medical Grade
- Warranty: Frame, Fabric, Padding - 2 years
- Regulatory: FDA listed, CE marked

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com